

Final Report for Medizinisches Labor Bremen

For Platelet Aggregation EQA Exercise Survey 2024-2

Exercise Date: Mittwoch, 9. Oktober 2024

Lab Code: 603

Report Generated: Mittwoch, 18. Dezember 2024 01:41

Number of Cases: 4

Number of Participants: 145

Report to Participants

We welcome members to participate then the platelet EQA survey 2024-2. The survey contains 4 platelet light transmission aggregation cases. From the history and platelet aggregation tracings, determine the following: 1. Agonist response based on reference ranges; all case platelet-rich plasma platelet counts were in the appropriate range or were adjusted unless otherwise noted. 2. Overall interpretation; choose the answer from the multiple choice list that most accurately states the interpretation. We also recommend that laboratories review the published North American guidelines for helpful hints on how to interpret aggregation findings (Hayward CM, Moffat KA, Raby A, et al: Development of North American Consensus for Medical Laboratories That Perform and Interpret Platelet Function Testing Using Light Transmission Aggregometry. Am J Clin Pathol 2010; 134:955-963)

